

# *The Best Just Got Better!*

## *Introducing the New* **EASYPULSE5+6** *Oxygen Conserving Regulator*



## FEATURES

### ***We've Added More Options***

The EasyPulse5+6 has 6 continuous flow settings, to compliment the 1 thru 5 pulse settings.

### ***Big Results in a Small Package***

Even with the added continuous flow settings, the EasyPulse5+6 weighs in at a mere 9.9 oz and is only 4.75 inches long, making it still the smallest and lightest conserving regulator available.



### ***Reduces Deliveries***

With the conserving ratio of up to 5.7 to 1, you enjoy the same great savings you have always enjoyed with the EasyPulse5.


### ***CMV Technology***


In the Pulse Mode, the EasyPulse5+6 delivers a front loaded tapered wave form. All the oxygen is delivered in the first third of the inspiratory phase. With the added security of more continuous flow options, you have the most versatile oxygen conserving device available.

**PRECISION**  **MEDICAL**

300 Held Drive • Northampton, PA 18067 USA • 800-272-7285  
[www.precisionmedical.com](http://www.precisionmedical.com)

**Made in**





## SPECIFICATIONS

### Dimensions:

Length 4.75" (12.07 cm)  
 Width 2.63" (6.67 cm)  
 Height 2.50" (6.35 cm)

### Weight:

9.9 oz (279 g)

### Inlet Pressure Range:

300-3000 psig

### Internal Regulated Pressure:

50 psig

### Altitude:

Sea level to 10,000 ft (3,048 m) above sea level

### Pulse Settings:

1, 2, 3, 4 and 5 L/min equivalents

### Flow Settings:

1, 2, 3, 4, 5 and 6 L/min Continuous

### Savings Ratio:

Up to 5.7 : 1

### Cannula Requirements:

Maximum 7 ft long adult standard single lumen oxygen nasal cannula

### Oxygen Cylinder Connection:

CGA 870 Post Valve

## ORDERING INFORMATION

Part#	Description
198712	EasyPulse5+6 Conserving Regulator with 1 Single Lumen Cannula (4 ft)
503943	M6 Cylinder with 870 Post Valve (empty)
504833	Single Lumen Cannula (4 ft)
504833-25	Single Lumen Cannula (case of 25)
503920	Carry Bag (fits M6 cylinder)
504184	Carry Bag (fits M4 and M6 cylinders)
504185	Carry Bag (fits ML6 and M9 cylinders)



ISO 13485  
 Certified

Precision Medical, Inc.  
 300 Held Drive • Northampton, PA 18067 USA  
 Phone: 610-262-6090 Toll-Free Phone: 1-800-272-7285  
 Fax: 610-262-6080 Toll-Free Fax: 1-800-352-1240  
 Visit us on the Web: [www.precisionmedical.com](http://www.precisionmedical.com)