

Mini Oxygen Regulator CGA 870 Style



"T" Handle
for easy
mounting
— no tools
needed

Barbed hose
outlet for easy
connection to the
oxygen tubing

Item# **RES042**

Adult • Barbed
0-8 LPM • **Each**
0, .5, 1, 1.5, 2, 2.5,
3, 4, 5, 6, 7, 8 LPM

Item# **RES043**

Adult • Barbed
0-15 LPM • **Each**
0, .5, 1, 1.5, 2, 3, 4,
5, 6, 8, 10, 15 LPM

Item# **RES046**

Pediatric • Barbed
0-4 LPM • **Each**
0, $\frac{1}{32}$, $\frac{1}{16}$, $\frac{1}{8}$,
 $\frac{1}{4}$, $\frac{1}{2}$, $\frac{3}{4}$, 1,
1 $\frac{1}{2}$, 2, 3, 4 LPM



HCPCS# E1353

For use with E size or smaller cylinders: A, B, C, D, JD, E, M2, M4, M6, ML6, M7, M9, M15, M22, or M24.